<p>Accu-CHEK® <i>Performa</i></p>
<i>Controls</i>
<p>REF 04861736</p>

®

Suitable for self-testing/For in vitro diagnostic use

Intended Use

For performance checks on the Accu-Chek Performa, Accu-Chek Performa Nano, Accu-Chek Performa Connect, Accu-Chek Performa Combo, Accu-Chek Performa Insight, and Accu-Chek Performa Solo systems with the Accu-Chek Performa test strips.

For performance checks on the Accu-Chek Performa (with code chip slot), Accu-Chek Performa Nano, and the Accu-Chek Inform II systems with the Accu-Chek Inform II test strips.

Testing control solutions with known glucose levels establishes that the operator and the system are performing acceptably. Control results must be within the defined acceptable ranges before valid patient testing is allowed.

Important Information

WARNING: Choking hazard. Small parts. Keep away from children under the age of 3 years.

DO NOT INGEST! Seek immediate medical attention if swallowed.

Contents of the pack

Pack containing 2 × 2.5 mL control solution and package inserts. The Accu-Chek Performa control solution pack contains two control solutions, one for the hypoglycaemic range (control solution 1, grey cap) and one for the hyperglycaemic range (control solution 2, white cap).

Disposal

All components of the pack can be discarded in domestic waste. Consult local ordinances as they may vary by country. Because the reactive substances are in such small quantities, they are not considered hazardous materials under EU regulations. If you have any questions, contact the local Roche representative.

Control solution storage and handling

- Refer to the meter User’s Manual for system operating conditions.
- The printed use by date is valid if the unopened control solutions are stored between 2 and 32 °C.
- The control solutions do not need to be kept in the refrigerator. DO NOT freeze the control solution.
- Write the date the control solution bottle was opened on the bottle label. The control solution must be discarded 3 months from the date the control solution bottle was opened (discard date) or on the use by date on the bottle label, whichever comes first.
- Control solutions taken directly from the refrigerator must be allowed to adjust to room temperature (without opening the control solution bottle).
- The control solution can stain fabric. Remove stains by washing with soap and water.

Performing a Control Test

Test control solutions the same way a blood sample is tested. See the meter User’s Manual for specific use instructions.

- Check the use by date on the test strip container. Do not use test strips past the use by date.
- Place the meter on a flat surface.
- Remove the cap from the control solution bottle. Wipe the tip of the bottle with a lint-free wipe.
- Squeeze the bottle until a tiny drop forms at the tip.
- Touch the drop to the **front edge** of the yellow window of the test strip.

The meter will indicate when there is enough control solution in the test strip.
- Wipe the tip of the bottle with a lint-free wipe. Cap the bottle tightly.

The control result appears on the display. Remove and discard the used test strip per facility policy.

You can compare the control result with the acceptable range printed on the test strip container label. If the control result is within the acceptable range, correct functioning of the system is assured.

If the control result is outside the acceptable range or if an error message is displayed, repeat the control test. If the second control result is also outside the acceptable range or an error message is displayed, contact the local Roche representative.

Sources of Error

If the control test results are outside the acceptable range, do not use the meter until you solve the problem. Check this list to help solve the problem.

- Have the test strips or control solutions expired?
- Did you wipe the tip of the control solution bottle before use?
- Were the caps on the test strip container and control solution bottle always closed tightly?

- Was the test strip used immediately after it was removed from the test strip container?
- Were the test strips and control solutions stored in a cool, dry place?
- Was the correct control solution level selected when the test was performed?
- Did you follow the directions?

Testing Intervals

Follow your facility’s policy for control testing intervals.

- Control tests should be performed:
- at the first time before using the meter for patient blood glucose testing.
 - at intervals established by the facility.
 - when a new test strip box is opened.
 - if the test strip container was left open.
 - if test strips were improperly stored.
 - if there is a question about a blood glucose result.
 - to check the performance of the system.
 - if the meter was dropped.

Your facility may require that control tests be successfully performed **after** any of the following occur and **before** patient testing resumes:

- Previous control results were outside the acceptable range.
- Control tests were not performed at the proper interval.

Control results must be within the designated range on the test strip container label, or as defined by your facility, before being considered acceptable. Patients can be tested after control tests have been acceptably performed at the proper testing interval.

Ingredients

Component	Control Solution Level 1 (Low)	Control Solution Level 2 (High)
Glucose	0.06 %	0.30 %
Buffer	4.83 %	4.82 %
Biological Salt	3.39 %	3.38 %
Preservative	0.30 %	0.30 %
Non-reactive ingredients	9.99 %	9.97 %
FD&C Blue #1	0.08 %	0.08 %
Percentages in weight/weight		

NOTE

For an explanation of symbols used refer to the end of the insert.

Visit our website at www.accu-chek.com or contact the local Roche representative for more information.

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